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Instructions for the reprocessing

(cleaning, disinfecting and sterilising)

of medical devices produced by Hager & Meisinger GmbH

Issued: January 2008

The medical devices produced and sold by Hager & Meisinger GmbH are re-usable unless their label contains explicit information to the contrary. However, as a rule, it is the sole responsibility of the doctor/expert using the devices to decide whether, depending on the respective case and the potential wear and tear of the products, he can re-use the products and how frequently he uses them. In case of doubt, it is always advisable to discard the products early and to replace them. The manufacturer Hager & Meisinger GmbH cannot guarantee the faultless function and performance of the products combined with a maximum degree of safety if the products are overused.

These reprocessing instructions apply in principle to all medical devices making up the product range supplied by Hager & Meisinger GmbH. Any particular features and/or exclusions that only concern individual items or groups of items are referred to separately. As to the general application and safety instructions concerning the use of the products, we would advise you to consult the application and safety instructions for the medical products of Meisinger, which are available separately (see also at www.meisinger.de).

General rules

All products need to be cleaned, disinfected and sterilised prior to use; this applies in particular to the first-time use of products i.e. after their delivery because all products are supplied non-sterile (cleaning and disinfecting after the removal of transport packaging; sterilisation after removing wrapping). The thorough cleaning and disinfecting is indispensable in order to achieve effective sterilisation.

As you are responsible for the sterility of the products during use, please ensure

- that, as a rule, only sufficiently validated device and product-specific procedures are used for the cleaning/disinfecting and sterilisation process,
- that the used devices (disinfector, sterilizer) are maintained and checked at regular intervals and
- that the validated parameters are adhered to during each cycle.

Please ensure that, during use, you separate the soiled products and do not return them to the bur block/instrument tray to avoid further contamination of the loaded bur block/instrument tray. After cleaning/disinfecting the soiled products, arrange them in the bur block/instrument tray and sterilise the fully loaded bur block/instrument tray. Please also take note of the statutory regulations valid in your country as well as the hygiene rules followed by the doctor's practice or the hospital. This applies in particular to the different rules regarding effective prion deactivation. You will find further instructions for the reprocessing of medical products and information on suitable products to carry out the reprocessing on the internet at www.dghm.de or at www.rki.de.

As the products are destined to be used for surgical, parodontological or endodontic procedures such as root canal debridement, they may penetrate the skin or the mucosa and come into contact with blood, internal tissues or organs (including wounds). Therefore, we recommend that they be assigned to risk group Critical B if used for their intended purpose.

Attention: For some products, additional aspects need to be considered (see chapter on 'Special instructions').

Cleaning and disinfecting

Basic rules:

Ideally, a mechanical procedure (disinfector) should be used for the cleaning and disinfecting process. Due to their clearly reduced effectiveness and reproducibility, manual procedures – even if they include the use of an ultrasound bath – should only be used if an automated procedure is unavailable. Pre-treatment needs to be carried out in both cases.

Apply hygiene-effective measures in line with the RKI guidelines at your workplace and carry out all work using powder and latex-free gloves.

Pre-treatment:

Abrasive impurities need to be removed from the products directly after use (within two hours maximum). To do so, use running water or a disinfectant solution; the disinfectant must not contain aldehydes (which could fix blood residues to the instrument surface), its effectiveness should be established (e.g. it should be certified by the VAH (German Association for Applied Hygiene), the DGHM (German Society for Hygiene and Microbiology) or the FDA or have CE marking), it should be suitable for the disinfection of the products and compatible with the products (see chapter on 'Material resistance'). For the manual removal of impurities, only use nylon brushes intended for the purpose (also see the chapter on 'Special Instructions') or clean soft and lint-free cloths that you only use for this purpose. Do not use metal brushes or steel wool.

As to products with lumen (cavities): rinse all cavities three times by using a disposable syringe (minimum volume 5-10 ml) and a cannula.

Please note that the disinfectants used during pre-treatment only ensure personal protection and can be no substitute for the disinfection procedure to be used later – after completion of the cleaning process.

Mechanical cleaning/disinfection (disinfector/RDG (cleaning and disinfection device):

- When choosing a disinfector you will have to ensure
- that the effectiveness of the disinfector has been certified (e.g. it has been licensed by the DGHM or the FDA or has CE marking according to DIN EN ISO 15883),
 - that, if at all possible, a programme is used that has been certified for thermic disinfection (A0-value > 3000 or – with regard to older devices – at least 5 minutes at 90 °C (194 °F)) (chemical disinfection runs the risk of disinfectant residues remaining on the instrument),
 - that the used programme is suitable for the products and has a sufficient number of rinsing cycles,
 - that only sterile water or water with low levels of germs (max. 10 germs/ml) and endotoxins (max. 0.25 endotoxin units/ml) is used for the post-purge cycle (e.g. purified water/highly purified water),

- that the air used for drying is filtered and that the disinfector is regularly maintained and checked.

When choosing an appropriate cleaning system you need to ensure:

- that it is generally suitable for the cleaning of products made of metal and plastic,
- that, in addition, – if no thermic disinfection is used – a disinfectant is used, whose efficacy has been certified (e.g. VAH/DGHM FDA licence or CE marking) and
- that it is compatible with the used cleaning agents and that the used chemicals are compatible with the products (see chapter on 'material resistance').

It is essential that the concentrations recommended by the manufacturer of the cleaning and the disinfectant agent (if required) are adhered to at all times.

Procedure:

1. Disassemble the products as far as possible.
2. When using products with lumen: ensure that all lumens are rinsed effectively as part of the pre-cleaning process.
3. Place the disassembled products in the disinfector. Take care that the products do not touch. If possible: connect all
4. Start the programme.
5. Remove the products from the disinfector after the programme has finished.
6. Check and wrap the products straight after removal if possible (see chapters on 'Checking', 'Maintenance' and 'Packaging') if necessary, after they have been dried off completely in a clean place.

Proof of the general suitability for effective mechanical cleaning and disinfecting has been provided by an independent certified test laboratory using the disinfector G 7836 GD (thermic disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the cleaning agent Neodisher medizym (Dr. Weigert GmbH & Co. KG, Hamburg). As part of the testing the laboratory used the above described procedure.

Manual cleaning and disinfection:

When choosing an appropriate cleaning and disinfecting agent you need to ensure

- that they are generally suitable for the cleaning and/or disinfection of products made of metal and plastic,
- that the cleaning agent, where used, is suitable for ultrasound cleaning (no production of foam)
- that an appropriate disinfectant is used whose efficacy has been certified (e.g. VAH/DGHM FDA licence or CE marking) and that it is compatible with the used cleaning agent and
- that the used chemicals are compatible with the products (see chapter on 'Material resistance').

Ideally, combined cleaning/disinfecting agents should not be used. Combined cleaning/disinfecting agents can only be used in cases where there is a very low degree of contamination (no visible soiling).

It is essential that the concentrations and contact times recommended by the manufacturer of the cleaning and the disinfectant agents are adhered to at all times.

Only use freshly made solutions and water that is either sterile or low in germs (max. 10 germs/ml) and endotoxins (max. 0.25 endotoxin units/ml) (e.g. purified water/highly purified water) and only use filtered air for drying.

Procedure: Cleaning

1. Disassemble the products as far as possible.
2. Place the disassembled products in the cleaning bath for the prescribed contact time and ensure that products are sufficiently submerged; where required, use ultrasound support but at least careful brushing using the nylon brushes intended for this purpose (see also chapter on 'Special Instructions'). Take care that the products do not touch. With regard to products with lumen: rinse all instrument lumens five times at the beginning and or at the end of the contact time using a disposal syringe (minimum volume 5-10 ml) and a cannula.
3. Afterwards, take the products out of the cleaning bath and rinse at least three times thoroughly with water. With regard to products with lumen: rinse all instrument lumens five times at the beginning and or at the end of the contact time using a disposal syringe (minimum volume 5-10ml) and a cannula.
4. Check the products (see chapters on 'Checking' and 'Maintenance').

Disinfection:

5. Place the disassembled, cleaned and checked products in the disinfecting bath for the prescribed contact time and ensure that the products are sufficiently submerged. Take care that the products do not touch. With regard to products with lumen: rinse all instrument lumens five times at the beginning and or at the end of the contact time using a disposal syringe (minimum volume 5-10 ml) and a cannula.
6. Afterwards, take the products out of the disinfecting bath and rinse them thoroughly with water at least three times. With regard to products with lumen: rinse all instrument lumens five times at the beginning and or at the end of the contact time using a disposal syringe (minimum volume 5-10 ml) and a cannula.
7. Dry the products by blowing them dry using filtered pressurised air.
8. Wrap the products, if possible, straight after removal (see chapter on 'Packaging', if necessary, after they have been dried off completely in a clean place).

Proof of the general suitability for effective manual cleaning and disinfecting has been provided by an independent certified test laboratory using the cleaning agent Gigazyme and the disinfectant Gigasept Instru AF (Schülke & Mayr GmbH, Norderstedt). As part of the testing, the laboratory used the above described procedure.

Checking

After all products have been cleaned and/or cleaned/disinfected, check them for corrosion, damaged surfaces/bare patches, broken/chipped-off edges, deformations (e.g. bent rather than round) and impurities and eliminate damaged products (limited numbers for re-use see chapter on 'Re-use'). Products that are still contaminated need to be cleaned and disinfected once more.

Maintenance

Re-assemble disassembled products (see specific instructions). Instrument oils must not be used.

Packaging

Arrange the cleaned and disinfected products in the dedicated bur block/sterilisation tray. Wrap the products and/or the bur blocks/sterilisation trays using disposable sterilisation packaging (disposable or double packaging) and/or sterilisation containers that meet the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilisation (thermal stability up to 137 °C (279 °F) sufficient steam permeability)
- Sufficient protection of the products and/or sterilisation packaging against mechanical damage
- Regular maintenance according to the manufacturer's specifications (sterilisation container)
- Individual packaging: the packaging must be sufficiently large to ensure that the sealing is tension-free.

Sterilisation

We only recommend the use of the sterilisation procedures listed below. Steam sterilisation

- Use of a fractional vacuum process or a gravitation process¹ (with sufficient product drying)
- Steam sterilizer according to DIN EN 13060 and/or DIN EN 285
- Validated according to DIN EN ISO 17665 (up to now: DIN EN 554/ANSI AAMI ISO 11134) (valid IQ/OQ (commissioning) and product-specific performance assessment (PQ))
- Maximum sterilisation temperature 134 °C (273 °F; plus tolerance according to DIN EN ISO 17665 (up to now: DIN EN 554/ANSI AAMI ISO 11134))
- Sterilisation time (exposure time at the sterilisation temperature) at least 20 min at 121 °C (250 °F) and/or at least 5² min at 132 °C (270 °F)/134 °C (273 °F)

¹ Use of the less effective gravitation procedure is only admissible if the fractional vacuum process is unavailable.

² or 18 min (prion deactivation)

Proof of the general suitability of the products for effective steam sterilisation has been provided by an independent certified test laboratory using the steam sterilizer Systec V-150 (Systec GmbH Labor-Systemtechnik, Wettengen) as well as the fractional process and the gravitation process. During the testing, the laboratory took into account the typical hospital and practice conditions as well as the procedure described above.

As a rule, the flash sterilisation procedure is not admissible. Furthermore, do not use hot air sterilisation, radiation sterilisation, formaldehyde or ethylene oxide sterilisation or plasma sterilisation. In order to avoid stains and corrosion, the steam must be substance-free (see limit values included in DIN EN 13060). When sterilising several devices, the maximum load of the sterilising apparatus must not be exceeded (observe the manufacturer's instructions).

Storage

Prior to the first use of the device, the product should be stored in its original packaging at room temperature in dust- and humidity-free conditions. Subsequently, the products should be stored in appropriate hygienically maintained containers (protected from dust, humidity and recontamination). After sterilisation, the products need to be stored in sterilisation wrapping in a dry and dust-free place. Please note the shelf-life resulting from the validation of the sterilisation wrapping.

Material resistance

When choosing the cleaning and disinfecting agents ensure that they do not contain the following components:

- organic, mineral and oxidising acids (minimum admissible pH value 5.5)
- strong alkaline solutions (maximum admissible pH value 8.5, neutral/enzymatic cleaning agent recommended)
- organic solutions (e.g. alcohols, ether, ketones, benzines)
- oxidants (e.g. hydrogen peroxides) - halogens (chlorine, iodine, bromine)
- aromatic/halogenated hydrocarbons

Never clean products, bur blocks and sterilisation trays using metal brushes or steel wool. Products, bur blocks and sterilisation trays should never be exposed to temperatures exceeding 137 °C (279 °F)!

Re-use

When it comes to products that are re-used, the end of their useful life tends to depend on their wear and tear as well as the use-related damage. Therefore, it is necessary to check the instrument after or prior to each re-use (see section on 'Checking'). Independently of this, it is the sole responsibility of the physician using the instrument to decide upon re-use based on the respective case as well as the potential wear and tear of the instrument. In this context, known restrictions as to the frequency of instrument use need to be observed.

If you fail to observe them, we do not accept any liability.

Special instructions

- Benex®-Control (Extractor):
- Disassemble instrument completely (including the Teflon plate) prior to cleaning and disinfecting
 - Reassemble instrument prior to wrapping and sterilisation; ensure that the ball does not catch in a cavity of the biteplate (semi-tilted position)
- Benex®-Control (Traction cable):
- Bend cable during pre-cleaning and manual cleaning and disinfection several times in all directions

Drill extension:

- Disassemble instrument completely prior to cleaning and disinfecting
 - Reassemble instrument prior to wrapping and sterilisation
- Diamond products and ceramic grinding tools:

- Use particular care during the cleaning of the grinding surfaces and ensure that all residues are removed

Distractors:

- Move instrument up and down several times during pre-cleaning and manual cleaning and disinfecting
 - Mechanical cleaning/disinfection as well as sterilisation only in specific positions (VCD05: half-open state, VCD15: nearly closed state)
- Products equipped with a conduit to supply a cooling medium and other products with lumens (canals, drillings etc.):
- Actively rinse the lumen during pre-cleaning and manual cleaning and disinfection process
 - Products containing lumens without feed-through channels must not be re-used

Supports:

- Disassemble instrument completely prior to cleaning and disinfecting

- Wrapping and sterilisation only in disassembled state.
- Bur blocks/instrument trays:
 - Cleaning and disinfecting only without products being loaded (products must not be cleaned and disinfected whilst they are in the bur block/instrument tray)
 - Remove attachments prior to cleaning and disinfecting and clean and disinfect in disassembled state
- Universal Torque Ratchet:
 - Disassemble instrument completely prior to cleaning and disinfecting
 - Move joint several times to and fro during pre-cleaning and manual cleaning and disinfecting
- Reassemble instrument prior to wrapping and sterilisation.

The above instructions have been validated by the manufacturer of the medical devices who found them to be SUITABLE for preparing a medical device for re-use. It is up to the person in charge of the reprocessing to ensure that, based on the use of the correct equipment, material and personnel in the reprocessing facility, the actual reconditioning process produces the desired results. Normally, this requires the validation and routine monitoring of the procedure. Equally, each deviation from the instructions provided should be carefully checked for effectiveness and potential adverse consequences by the person in charge of reprocessing.

General application and safety instructions for the medical device produced by MEISINGER

Issued: January 2008

- MEISINGER medical products (dental, implantology, Bone Management®, maxillary surgery, general surgery, chiropody devices) must only be used by dentists, doctors and/or the respective experts who, due to their training and experience, are intensely familiar with the use of these products and who have the corresponding expertise in the respective specialist fields. The use of surgical products requires relevant expertise and experience in dental implantology, maxillary surgery and/or other surgical fields including diagnosis, pre-operative planning and surgical techniques.
- It is the sole responsibility of the doctor in charge who, depending on the respective situation (indication), decides on the actual use of the MEISINGER products for each individual case.
- Courses providing further training on the possible uses of MEISINGER products are organised on a regular basis. You will find relevant information about these events as well as further information on Meisinger products on the internet at www.meisinger.de.
- All Meisinger products have been developed for specific applications. Therefore, inappropriate use can result in the premature wear and tear of the products and put patients and users at risk.

Application

- It is essential to only use turbines as well as hand and angle pieces that are technically and hygienically faultless, maintained and cleaned.
- Rotating instruments need to be clamped as far down as possible with their speed set before applying them on the object. are used with the rotary instruments.
- Using the instruments for canting or leveraging should be avoided as it increases the risk of breakage.
- Depending on the application, it is recommended to use protective goggles while using the instruments. Users of diamond disks should use a disk protection device.
- Inappropriate use of the products leads to badly executed work and increased risk.
- When working with dry materials, it is recommended to use a suction cleaning device.
- In particular, users of hand tools should take care to use them gently and with consideration.
- The user must at all times avoid touching the instruments and parts without protection (protective gloves should be worn).
- Thermal bone damage caused by rotating tools (e.g. pilot burr; conical burr, expansion burr) must at all times be avoided (user training, working at low speed and with sufficient cooling).

Use of pressure

- Users of the instruments should at all times avoid applying excessive pressure. This can damage the working part of the instruments and cause the cutting edges to break off. At the same time, it generates excessive heat.
- The use of excessive pressure when using grinding tools can cause the abrasive particles to break off or the instrument to become clogged and lead to heat generation.
- During polishing, excess pressure can lead to heat generation.
- Due to overheating, excess pressure can damage the dental pulp or, due to broken off cutting edges, it can result in undesired rough surfaces. In such cases, even instrument breakage cannot be excluded.

Cooling

- In order to avoid excessive heat generation during preparation, a sterile water/sodium chloride solution supplied via a permanent external feeding device should be used to ensure sufficient cooling during use of the instruments.
- When using FG instruments that are more than 22 mm long or whose head diameter exceeds 2 mm, additional external cooling is required.
- Insufficient cooling will lead to irreversible damage to the bone and/or the adjacent tissue.

Storage, disinfection, cleaning and sterilization

- Unless there is explicit information to the contrary, all MEISINGER products are supplied in non-sterile packaging and, depending on the application, they need to be sterilised prior to use. Prior to their first use on the patient and immediately after each use, all products need to be disinfected and sterilised. Inappropriate cleaning and sterilising of the instruments can result in the patient being infected with harmful bacteria.
- You will find detailed instructions for the disinfecting, cleaning and sterilising of products in the Instructions for the reprocessing of medical devices produced by Hager & Meisinger GmbH. We would be happy to provide you with these instructions at your request. They are also available on the internet at www.meisinger.de.
- The products should be stored in appropriate, hygienically maintained containers. The same applies to sterilised instruments. The stored products must be protected from dust, humidity and recontamination. Instructions as to maximum storage duration must be adhered to.

Speed recommendations

- Following the instrument-specific speed recommendations produces the best results.
- Exceeding the maximum admissible speed (rpm) when using long and pointed instruments tends to produce vibrations that can lead to the destruction of the instrument.
- When using working parts with diameters exceeding the thickness of the shaft, excessive speed can release great centrifugal forces that may cause the shaft to bend and/or the instrument to break. Therefore, the maximum admissible rpm must never be exceeded.
- Please consult the manufacturer's information (see catalogue or

www.meisinger.de) for the recommended and the maximum admissible speed ranges. Non-compliance with the maximum admissible speed puts safety at risk.

- Generally, the following rules apply:
 - The larger the working part of an instrument the lower the speed
 - The larger the working part of an instrument the greater the force of pressure
- Instruments with a maximum speed of 300,000 rpm are suitable for micro-motor hand pieces and turbines with stable ball bearings. Not to be recommended for turbines with air bearings.
- Instruments with a maximum speed of 30,000 to 160,000 rpm are suitable for micro-motor hand pieces or technical hand pieces up to the indicated rpm. Not to be recommended for turbines.
- Surgical instruments: suitable for geared down micro-motor hand and angle pieces 1:10 with stable ball bearings. Speed 600 to 800 rpm with physical and, possibly, sterile external cooling or internal cooling when using the respective hand piece.

Discarding worn instruments and parts

- As a rule, MEISINGER medical products are re-usable unless there is explicit proof and information to the contrary. It is the sole responsibility of the doctor using the devices to decide on the possibility of and the justification for multiple use as well as on the frequency of use, based on the respective case and the potential wear and tear of the products. If in doubt, it is always best to discard the products early and to replace them.
 - Instruments are subject to wear and tear. This is particularly true for rotating tools. Therefore, our warranty for all rotating tools expires after their fifth use. If the tools are used more frequently, the manufacturer, Hager & Meisinger GmbH, can no longer guarantee a faultless cutting performance combined with minimum temperature change.
 - Broken off cutting edges of instruments cause vibrations and great forces of pressure, which, in turn, leads to broken preparation corners and rough surfaces.
 - Bare patches on diamond instruments indicate a lack of abrasive particles and can be a sign of blunt instruments. This leads to excessive temperatures during instrument use.
 - Instruments that are bent and/or do not run true should be discarded forthwith.
-  • With the reuse of disposable products the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed

Additional instructions for the use of implants

- As a rule, implants are meant to be used only once.
- Implants must only be inserted by dentists and doctors who are intensely familiar with dental implantology including diagnosis, pre-operative planning, surgical techniques and prosthetic care.
- The unprotected touching of the implant by the user must be avoided (protective gloves must be worn).
- Systemic errors that are due to a mix-up of tools and implants are to be avoided. Therefore, paying attention to colour coding and labelling is essential.
- A lack of primary stability can be due to the faulty preparation of the implant bed as well as to deficient bone quality. In such cases, the mobile implant should be removed and replaced with a larger one.
- The best way to avoid complications is to only use implants when all the conditions for a complication free treatment guaranteeing long-term success have been met. If, from the start, there are reservations regarding the chances of success of a desired implant-based Dental prosthesis, it may be best to choose another prosthetic solution that does not require an implant.

Additional instructions regarding the use of trepans

- When using trepans, you have to proceed with particular care. For example, it is advisable not to exceed the recommended rpm speed ranges.
- In order to prepare for the actual use of a trepan, it should be set to produce counter-clockwise rotations creating a groove in the bone. Afterwards the trepan can be inserted into this groove and, using clockwise rotations, it can be moved further down.
- Carrying out a prior X-Ray is essential to establish the maximum possible drilling depth and to maintain the necessary distance, for example, to the mandibular nerve. As an additional safety measure to spare the nerve, the axial direction of the trepan countersink attachment, based on the sagittal level of the ascending branch, must be milled laterally at an angle of approx. 15-20°.

Further comments:

- Due to statutory regulations, returned goods can, on principle, only be accepted if the complete batch number is provided. This number can be found on the product packaging.