STEADY- RESIN – Information

1. Description

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STEADY-RESIN – Product Information

1. Description

Steady-Resin is a highly transparent cold cure resin, especially developed for the orthodontic field and is available in a polymer type S, M and VARIO – manufactured for different procedures (pls. see below). It is fabricated on the basis of methyl-methacrylate, is free of cadmium and contains tertiary amines for polymerization. When used properly, the biological tolerance is in accordance with the guide lines EN 30993 and EN/ISO 7405 (a certificate is available on request). It is exclusively designed for dental field application: e. g. expansion plates, retainers, activators, splints with articulated surfaces, etc. Colour dyes are available for the monomer as well as colour glitters for the polymer. Caution! Too much dye and glitter can effect the polymerisation resp. fit.

2. Hints for Use / Application

2.1 Preparation

It is recommended to soak the models in a 35 °C – 40 °C warm water for approximately 10 minutes in order to avoid bubbles on the basal surfaces.

2.2 Insulating

Coat the model once with a thin layer of orthodontic insulating agent (REF 8364).

2.3 Indication

Steady-Resin S – for the “salt and pepper” technique – has an excellent stability. It is applied directly onto the model which is insulated by the orthodontic insulating agent; starting with monomer, polymer and monomer are applied alternately by using spray bottles (REF 8141) and spray nozzles (REF 8138/8139). The finishing layer should be monomer again.

Steady-Resin M – for the “premix” technique – has a top pouring quality and a broad range of working possibilities. It can be mixed in a Resimix®-cup and after that applied on the insulated model.

Steady-Resin VARIO – Due to an extended setting time, Steady-Resin VARIO can be used for both, the “salt and pepper” and “premix” technique.

2.4 Working Parameters

The recommended times and temperatures are approximate and can be affected by room temperature and mixing ratio.

Handling Time: 6 – 8 min.
Mixing Ratio:
Steady-Resin »S« 10 : 4 weight percent
(polymer : monomer)
Steady-Resin »M«/
Steady-Resin VARIO 10 : 5 weight percent
(polymer : monomer)
This mixing ratio is approximately
2,7 : 1 volume percent
polymer : monomer.

2.5 Polymerization

Curing is strictly effected in a pressure pot filled with clean water. The procedure should not be interrupted and keep the values mentioned below. Do not exceed or fall below a value. Changes of parameter can lead to bubbles, inhomogeneity, heavy shrinkage and a higher content of monomer.
Temperature: approx. 45 - 48 °C / ≈ 113 °F
Pressure: 1.8 - 2.0 bar / ≈ 30 p.s.i.
Time: 15 - 20 min.

3. Finishing

Steady-Resin can be finished with carbide burs after the polymerization and polished conventionally. As with all methyl-methacrylate products it is recommended to provide a suction unit. For minimizing the remaining parts of monomer it is recommended to place it in water for several hours for hypersensitive patients.

4. Durability / Storage

Polymer is durable for 3 years at least in original containers. Monomer can be stored in original containers and for at least 24 months. It can be used after the durability when being checked previously. After pouring into transparent dosing bottles, keep them out of light as otherwise reactions on polymerization and changes in colour may occur. Store closed containers only in a cool (max. 25 °C) and dry place! Keep out of light!

5. Danger Hints

Liquid is highly flammable! Causes irritations on eyes, skin and respiratory organs. Keep away from sources of ignition and store only in well-ventilated rooms, do not allow to enter into drains. Avoid longer skin contact with non-polymerized material and inhaling of monomer vapours. In individual cases allergenic reactions on components of Steady-Resin can occur (e. g. methyl-methacrylate, N-Dimethyl-p-toluidine or dibenzoylperoxid). Furthermore, please pay attention to our Safety Data Sheet for Steady-Resin monomer.

All hints for the use of our materials in oral, written or practical form are given to the best of one's knowledge and shall be understood as hints. The use happens beyond our control and is subject to the responsibility of the user.

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